

Non-invasive Online CO₂ Measurement with Chemical optical CO₂ Sensors

Chemical optical CO₂ sensors provide a versatile tool for online monitoring of dissolved CO₂ in various scale and designs

Athanas Apostolidis, PreSens Precision Sensing, 93053 Regensburg, Germany

Accurate CO₂ measurement in the liquid phase has been a challenge for decades. Although off-gas measurement by means of IR-spectroscopy is a proven technology, the correlation of CO₂-levels in gas and liquid phase in biological systems is difficult.



Fig. 1 pCO₂ mini fiber optic CO₂ transmitter

The pCO₂ mini transmitter is a precise fiber optic CO₂ transmitter. It is used with CO₂ sensors based on a 2 mm optical fiber. A PC is connected to run the easy-to-use software. pCO₂ mini is additionally equipped with an analogue out and a trigger input. The device is also equipped with a PT 1000 temperature sensor for online temperature measurement and compensation of the measured signal.




Fig. 2 pCO₂ view v1 measurement control software for data acquisition, display and storage

The pCO₂ view software is used for controlling the operation of the pCO₂ mini transmitter. The data sent from

the device are shown in digital display and graphically in a trace chart. The software includes menus for customer specific calibration of sensors as well as the possibility to operate with predefined calibration values. Data can be stored any time within a measurement and exported to MS Excel or plain ASCII-files. The upload of stored measurements is also possible for later evaluation.

CO₂ sensor performance

The CO₂ sensors are optimized for physiological solutions. They are calibrated with gas mixtures of defined CO₂ concentration levels and delivered with a final inspection protocol. The customer can use either the factory data for calibration of the pCO₂ mini transmitter or perform his / her own calibration. The dynamic range of the CO₂ sensor is 1 % to 25 % CO₂ in order to fit to typical CO₂ levels within cell culture or other biotechnology applications. Sensors can be delivered beta-irradiated or untreated. Autoclavable sensors are also available.




Fig. 3 Time trace of sensor reading at 20 °C upon exposure to atmospheres with defined levels of up to 30 % CO₂ in the mixture.

Measurement principle

The chemical optical carbon dioxide sensor is based on our patented DLR measurement principle^[1]. The light of the blue LED excites the sensor to emit luminescence. The luminescence lifetime measured is a superposition of the signals of an analyte sensitive indicator and an inert reference indicator, where both indicators exhibit very

CO₂ Measurement

different luminescence lifetimes and the luminescence of the CO₂ sensitive indicator can be suppressed by CO₂ (see [2] in Fig. 4). The measurement signal correlates to the partial pressure of carbon dioxide. This internal referencing compensates intensity fluctuation caused e. g. by changes in excitation intensity.

Bioreactor Monitoring and Control. ESACT 2009, Dublin, Ireland




Fig. 4 Measurement principle of CO₂ sensors: [1] signal in absence of CO₂, [2] signal of sensor in CO₂ containing environment

PreSens CO₂ sensors as an alternative to off-line sampling [2]

PreSens CO₂ sensors were evaluated for monitoring the cultivation cells in a 3-L bioreactor. As shown in Fig. 5, the sensors show excellent long-term stability over a period of 10 days during cell cultivation of CHO cells with an initial seed density of 0.4×10^6 cells/ml. The readings of the PreSens CO₂ sensor were also in good correlation to off-line CO₂ measurements. Additionally, the PreSens CO₂ sensor offers the possibility to continuously monitor the CO₂ concentration between two off-line measurements.




Fig. 5 Comparison of on-line (PreSens CO₂ sensor) and off-line measurement of CO₂ in a 3-L bioreactor

References:

- [1] Patent DLR Nr. EP 1000345 [B1] and Nr. US 6602716 [B1]
- [2] Baggio, et al. Physical-Chemical Sensors: Application to

Bring to light what's inside. Ask our experts:

PreSens Precision Sensing GmbH
 Josef-Engert-Str. 11
 93053 Regensburg, Germany

Phone +49 941 94272100
 Fax +49 941 94272111
 info@PreSens.de

www.PreSens.de