

# Real Time Respiration Measurements

## Measuring the aerobic respiration of complex microbial communities

B. R. Roller<sup>1</sup>, Z. M. Lee<sup>1</sup>, K. E. Studer-Rabeler<sup>2</sup>, and T. M. Schmidt<sup>1</sup>

<sup>1</sup>Dept. of Microbiology and Molecular Genetics, Michigan State University, East Lansing, MI, USA

<sup>2</sup>Coy Laboratory Products, Grass Lake, MI, USA

**Many methods for measuring aerobic respiration, including electrodes and colorimetric assays, limit the capacity for obtaining replicate measurements or the ability to detect oxygen consumption in real time. The OxoDish® - a microtiter plate with integrated oxygen sensors - and SDR SensorDish® Reader offer the potential for making twenty four simultaneous respiration measurements in real time, benefitting ecologists measuring the respiration of complex microbial communities.**

Our aim was to optimize the SDR OxoDish® system for measuring the rate of oxygen (O<sub>2</sub>) consumption of pure cultures of bacteria and in soil slurries of complex microbial communities. Soil slurries are particularly challenging due to the potential for physical settling of soil particles on the optode. Two modifications of the OxoDish® were needed to measure O<sub>2</sub> consumption accurately for these samples: minimizing O<sub>2</sub> influx into the wells of the microtiter plate and maintaining a well-mixed solution so that O<sub>2</sub> consumption was homogeneous throughout the depth of the well. A pure culture of *Pseudomonas* strain HF3 and forest soil were used to test modifications to the OxoDish®.

### Materials & Methods

The influx of O<sub>2</sub> into wells covered with adhesive foil was determined experimentally. A growth medium containing Resazurin, a colorimetric indicator of O<sub>2</sub>, and cysteine, a reducing agent that consumes O<sub>2</sub>, was added into 12 wells of an OxoDish® under anaerobic conditions in a vinyl anaerobic chamber / glove box. The wells were filled as close as possible to capacity and half of the wells were covered with adhesive foil and compared to uncovered wells in an ambient oxygen atmosphere at 25 °C and 100 RPM on a rotary shaker. The effect of glass beads, used to

promote mixing beyond rotary shaking, on O<sub>2</sub> consumption was tested with sterile growth medium as the negative control. Samples of *Pseudomonas* strain HF3 were removed from a balanced growth culture and placed into six OxoDish® wells, half of which contained 2 sterilized glass beads per well, and covered with adhesive foil (Fig. 1). O<sub>2</sub> concentration was measured for each well every 15 seconds for 120-minute experiments at 25 °C and using rotary shaking at 100 RPM. In all rate measurements, the first 5 minutes of collected data were trimmed to permit time for equilibration.

The bead and foil modifications to the OxoDish® were then utilized with soil slurry samples. Forest soil was sieved (2 mm), and resuspended 1 : 1 (w : v) in 10 mM MES buffer (pH 6) supplemented with 50 µg/mL cycloheximide and 2.3 mM sodium pyrophosphate. The soil slurry was stirred for 30 minutes, filtered through 8 layers of cheese cloth and strained through either a 100 µm or 40 µm cell strainer for different particle loading. The slurry was then added to wells with two sterile beads as described for pure culture experiments and sterile buffer was used as negative control.

### Oxygen Influx

Wells that contained sterile anaerobic medium and were covered with adhesive foil showed no detectable O<sub>2</sub> influx for more than 120 minutes, while uncovered wells began showing O<sub>2</sub> influx after 40 minutes of exposure to an ambient atmosphere while shaking at 100 RPM (Fig. 2). The 40 minute delay in SDR detection of O<sub>2</sub> influx of uncovered wells can be attributed to cysteine consuming any O<sub>2</sub> that diffused into the media. Resazurin confirmed the SDR readings, with 6 uncovered wells turning pink within minutes of removal from the anaerobic environment, indicating O<sub>2</sub> influx, while 5 of the 6 covered wells remained clear, indicating no O<sub>2</sub> influx over many hours.




Fig. 1: Image of the OxoDish® with bead and foil modifications.

# SDR Respiration




Fig. 2: The effectiveness of adhesive foil at limiting  $O_2$  influx into wells containing anaerobic growth medium.




Fig. 4: The effect of decreased particle loading on  $O_2$  consumption rates.

## Pure Culture Oxygen Consumption

The second experiment measured the respiration of *Pseudomonas* strain HF3 and the effect of mixing with beads. Bacterial  $O_2$  consumption in wells that were mixed with beads was less variable than wells without beads, which is attributed to a more homogeneous solution throughout the well depth (Fig. 3). Using linear regression, the mean  $O_2$  consumption rate for samples with beads was 5.74 nmol/minute and replicates had a coefficient of variation (c. v.) of 6.7%. The mean rate for samples without beads was 11.86 nmol/minute (c. v. = 18.9%). The rates of  $O_2$  consumption have not been corrected for ingress of  $O_2$ , but they are adequate for comparisons of variability between treatments. There was no detectable effect of cells settling on the optode when stirred with beads, as the rate of  $O_2$  consumption doubled after one doubling time of the pure culture population (data not shown). This indicates the solution was homogeneous and not experiencing a shift in growth when transferred from a flask to a sealed OxoDish® well. The specific  $O_2$  consumption rate with beads (mean  $2.0 \times 10^{-7}$ ), while not corrected for ingress, is similar to the rate measured by microrespirometry (mean  $3.4 \times 10^{-7}$ ).




Fig. 3: *Pseudomonas* strain HF3  $O_2$  consumption during balanced growth.

## Soil Slurry Oxygen Consumption

The beads and foil modifications to the OxoDish® were then extended for use with a complex microbial community.  $O_2$  consumption rates were measured for a soil slurry strained through two different sized cell strainers for different particle loading. The biggest difference between the two particle loading treatments came in the variability of measurement, with the 40  $\mu\text{m}$ -strained samples having a coefficient of variation for the mean  $O_2$  consumption rate of 7.08%, which is much smaller than the 100  $\mu\text{m}$ -strained samples coefficient of variation of 35.05% (Fig. 4). The decrease in variability for 40  $\mu\text{m}$ -strained samples is attributed to increased homogenization of the soil slurry in the well, with the 100  $\mu\text{m}$ -strained samples having too many particles that settled.

## Conclusion

The SDR used with an OxoDish® is useful for the microbial ecology research community due to its capacity for replication and real time measurements. The OxoDish® can be modified with adhesive foil to minimize  $O_2$  influx into the wells of a microtiter plate and with glass beads to maintain homogeneous solutions throughout the depth of the well. Future investigations of complex microbial communities using the SDR and OxoDish® can measure  $O_2$  consumption while manipulating multiple variables simultaneously, where a large capacity for replication is necessary.

Bring to light what's inside. Ask our experts:

PreSens Precision Sensing GmbH  
Josef-Engert-Str. 11  
93053 Regensburg, Germany

Phone +49 941 94272100  
Fax +49 941 94272111  
info@PreSens.de

[www.PreSens.de](http://www.PreSens.de)