

Scientific Paper:

Biotechnology Letters 26, 763-767, 2004

On-line oxygen uptake rate and culture viability measurement of animal cell culture using microplates with integrated oxygen sensors

Rahul Ravi Deshpande & Elmar Heinzle*

Biochemical Engineering Institute, Saarland University, 66123 Saarbruecken, Germany

*Author for correspondence (Fax: +49-681-3024572; E-mail: e.heinzle@mx.uni-saarland.de)

Abstract:

O₂ uptake rates of animal cells (Chinese hamster ovary-CHO) were measured in 96-well microtiter plates by integrating with fluorescent sensors thereby measuring fluorescence intensity ratios of an O₂-sensitive and an insensitive fluorophor. O₂ consumption rate was estimated from measured dissolved O₂ and from O₂ mass transfer coefficient determined in advance. Specific uptake decreased with time from 3.2×10^{-13} mol O₂ cell⁻¹ h⁻¹ at 15 h cultivation to 1.8×10^{-13} mol O₂ cell⁻¹ h⁻¹ at 48 h. Specific O₂ uptake was also determined by sampling from a spinner-flask culture giving identical values. A cell viability assay for cultures based on O₂ measurements is described in which cells are incubated outside the fluorescence reader and then the dissolved O₂ is measured only once at a fixed time after the start of incubation. This protocol can be directly applied for high-throughput measurements.

Key-words: CHO cells, culture viability, fluorescent oxygen sensor, microtiter plates, oxygen uptake rate