

**EOM-pH-mini****Specifications**

pH sensor	HP5 / HP8
Phase resolution	0.05 °
Minimal sampling rate	Below 1 sec.
Supply voltage	7 - 18 V
Power consumption in active mode	1.5 W
Power consumption in stand-by mode	0.6 W
Temperature: operating / storage	0 °C to 50 °C / - 10 °C to 70 °C
Relative humidity	Up to 80 % (non condensing)
Dimensions	40 mm x 100 mm x 22 mm
Weight	70 g
Digital interface	RS232 interface (optional RS485)