

| | OXY-4 micro | OXY-10 micro |
|----------------------------------|--|--------------|
| Specifications | | |
| Oxygen sensors | PSt1 only (optical ST connector) | |
| Power supply | 110 - 230 VAC, 50/60 Hz, max. 40 vA (cables included) | |
| Temperature: operating / storage | 0 °C to 50 °C / - 10 °C to 60 °C | |
| Relative humidity | Up to 80 % (non condensing) | |
| Dimensions | 300 mm x 250 mm x 135 mm | |
| Weight | 3.3 kg | 3.9 kg |
| Digital interface | RS232 interface (with RJ connector to serial port, cable included) | |