

	OXY-4 mini & OXY-10 mini	OXY-4 trace & OXY-10 trace
Specifications		
Oxygen sensors	PSt3 only (optical SMA connector)	PSt3 and PSt6 (optical SMA connector)
Power supply	110 - 230 VAC, 50/60 Hz, max. 40 vA (cables included)	
Temperature: operating / storage	0 °C to 50 °C / - 10 °C to 60 °C	
Relative humidity	Up to 80 % (non condensing)	
Dimensions	300 mm x 250 mm x 135 mm	
Weight OXY-4 mini & OXY-4 trace	3.3 kg	
Weight OXY-10 mini & OXY-10 trace	3.9 kg	
Digital interface	RS232 interface (with RJ connector to serial port, cable included)	