

	LP-1	Battery Cable	Single USB Serial Adapter/ 4-Port USB Serial Adapter
<b>Specifications</b>			
Power Supply	Provided by connected instrument (e.g. Fibox 3), Input voltage range: 9 – 24 VDC, 5 W (nominal 18 VDC)	Input voltage range: 9 – 24 VDC, 5 W (nominal 18 VDC)	No external power adapter required, powered by USB
Compatibility	Microx TX3 (trace) Fibox 3 (trace v3)	o For single-channel transmitters: Microx TX3 (trace), Fibox 3 (trace), pH-1 mini, pH-1 micro o For LCD transmitters: Fibox LCD (trace)	Windows 98/ME/2000/XP/Vista (Certified Microsoft WHQL approved), driver CD include
Dimensions (D x W x H) / Weight	190 x 125 x 45 mm / 550 g	Length: 2.0 m / 100 g	USB cable length: 1.0 m (included), up to 5 m detachable
Oxygen Units	% air saturation, % O <sub>2</sub> , hPa, Torr, mg/L (ppm), μmol/L		
Data Logger	More than 20,000 samples stored in internal EEPROM		
Graphic Display	128 x 64 pixels, backlight, automatic contrast control		
Operating / Storage Temperature	0 °C to 50 °C / - 10 °C to 65 °C		
Relative Humidity	Up to 60 % (IP20)		

	Adapter for Round Containers (ARC)	Stick-On Adapter (SOA)
<b>Specifications</b>		
Dimensions (D x W x H)	Velcro <sup>®</sup> strip 1000.0 x 22.0 x 4.0 mm	20.0 x 20.0 x 7.0 mm
Connector Type	SMA socket	SMA socket
Compatible with All Devices With SMA Connectors	Fibox 3, pH-1 mini, OXY-4 mini, OXY-10 mini, and others	Fibox 3, pH-1 mini, OXY-4 mini, OXY-10 mini, and others

	Coaster for Shake Flasks (CFG) / Centrifuge Tube Adapter (CTA)	POF
<b>Specifications</b>		
Dimensions (D x W x H)	Approx. 93.0 x 41.0 x 16.0 mm	Optical diameter is 2 mm; outer diameter including the black cladding is approx. 2.8 mm
Connector Type	SMA socket	SMA connectors on one or both sides available for use with stick-on adapter and adapter for round containers
Length of Fiber	2.5 m	Available lengths for the POF are 1.0, 2.5 and 5.0 m (for lengths of more than 5 m, please contact our service team)
Compatibility	All devices with SMA connectors, e.g. Fibox 3, pH-1 mini, OXY-4 mini, OXY-10 mini	All devices with SMA connectors, e.g. Fibox 3, pH-1 mini, OXY-4 mini, OXY-10 mini
Details	Compatible with shake / spinner flasks up to 1l	Temperature stability: The POF is resistant to temperatures up to 70 °C