

Specifications	OXY-4 micro	OXY-10 micro
Oxygen sensors	PSt1 only (optical ST connector)	
Power supply	110 - 230 VAC, 50/60 Hz, max. 40vA (cables included)	
Temperature: operating / storage	0°C to 50 °C / -10°C to 60°C	
Relative humidity	up to 80 % (non condensing)	
Dimensions	300 mm x 250 mm x 135 mm	
Weight	3.3 kg	3.9 kg
Digital interface	RS 232 interface (with RJ connector to serial port, cable included)	